

Amaro 5000 ORIZONTAL

Modele	Dimensiuni interioare	Volum	Dimensiuni exterioare	Unitati de sterilizare
	Cm	Litri	Cm	STU
70.112.80	70x112x80	622	208x221x120	6
70.112.110	70x112x110	856	208x221x150	9
70.112.140	70x112x140	1090	208x221x180	12
70.112.170	70x112x170	1323	208x221x210	15
70.112.200	70x112x200	1557	208x221x240	18
70.112.230	70x112x230	1790	208x221x270	21
70.180.130	70x180x130	1638	208x256x166	20
70.180.160	70x180x160	2016	208x256x196	25
70.180.190	70x180x190	2394	208x256x226	30
70.180.250	70x180x250	3150	208x256x286	40

